

## Lab 4: Mohr Titration of Sodium in Pickles

### Questions

1. Standard A:  $\frac{0.1034 \text{ g NaCl}}{42.42 \text{ mL} - 24.30 \text{ mL}} \times \frac{1000 \text{ mL}}{1 \text{ L}} \times \frac{1 \text{ mol NaCl}}{58.443 \text{ g NaCl}} \times \frac{1 \text{ mol AgNO}_3}{1 \text{ mol NaCl}} = 0.09764 \text{ M AgNO}_3$

Standard B:  $\frac{0.0992 \text{ g NaCl}}{26.15 \text{ mL} - 9.18 \text{ mL}} \times \frac{1000 \text{ mL}}{1 \text{ L}} \times \frac{1 \text{ mol NaCl}}{58.443 \text{ g NaCl}} \times \frac{1 \text{ mol AgNO}_3}{1 \text{ mol NaCl}} = 0.100 \text{ M AgNO}_3$

Standard C:  $\frac{0.0982 \text{ g NaCl}}{43.52 \text{ mL} - 26.39 \text{ mL}} \times \frac{1000 \text{ mL}}{1 \text{ L}} \times \frac{1 \text{ mol NaCl}}{58.443 \text{ g NaCl}} \times \frac{1 \text{ mol AgNO}_3}{1 \text{ mol NaCl}} = 0.0981 \text{ M AgNO}_3$

Average:  $\bar{x} = \frac{\sum x_i}{N} = \frac{0.09764 \text{ M AgNO}_3 + 0.100 \text{ M AgNO}_3 + 0.0981 \text{ M AgNO}_3}{3} = 0.099 \text{ M AgNO}_3$

Standard Deviation:  $s = \sqrt{\frac{\sum (x_i - \bar{x})^2}{N-1}} =$

$$\sqrt{\frac{(0.09764 \text{ M AgNO}_3 - 0.099 \text{ M AgNO}_3)^2 + (0.100 \text{ M AgNO}_3 - 0.099 \text{ M AgNO}_3)^2 + (0.0981 \text{ M AgNO}_3 - 0.099 \text{ M AgNO}_3)^2}{3-1}} = 0.001 \text{ M AgNO}_3$$

2. Pickle Aliquot 1:  $(30.37 \text{ mL} - 11.20 \text{ mL}) \times 0.099 \text{ M AgNO}_3 \times \frac{1 \text{ L}}{1000 \text{ mL}} \times \frac{1 \text{ mol Na}}{1 \text{ mol AgNO}_3} \times \frac{22.990 \text{ g Na}}{1 \text{ mol Na}} \times \frac{1000 \text{ mg}}{1 \text{ g}} = 43. \text{ mg Na}$

Pickle Aliquot 2:  $(20.58 \text{ mL} - 0.98 \text{ mL}) \times 0.099 \text{ M AgNO}_3 \times \frac{1 \text{ L}}{1000 \text{ mL}} \times \frac{1 \text{ mol Na}}{1 \text{ mol AgNO}_3} \times \frac{22.990 \text{ g Na}}{1 \text{ mol Na}} \times \frac{1000 \text{ mg}}{1 \text{ g}} = 44. \text{ mg Na}$

Pickle Aliquot 3:  $(40.27 \text{ mL} - 20.84 \text{ mL}) \times 0.099 \text{ M AgNO}_3 \times \frac{1 \text{ L}}{1000 \text{ mL}} \times \frac{1 \text{ mol Na}}{1 \text{ mol AgNO}_3} \times \frac{22.990 \text{ g Na}}{1 \text{ mol Na}} \times \frac{1000 \text{ mg}}{1 \text{ g}} = 44. \text{ mg Na}$

Average:  $\bar{x} = \frac{\sum x_i}{N} = \frac{43. \text{ g Na} + 44. \text{ g Na} + 44. \text{ g Na}}{3} = 44 \text{ g Na}$

Standard Deviation:  $s = \sqrt{\frac{\sum (x_i - \bar{x})^2}{N-1}} =$

$$\sqrt{\frac{(43. \text{ g Na} - 44 \text{ g Na})^2 + (44. \text{ g Na} - 44 \text{ g Na})^2 + (44. \text{ g Na} - 44 \text{ g Na})^2}{3-1}} = 1 \text{ g Na}$$

3.  $\frac{43. \text{ mg Na}}{20.00 \text{ mL pickle water}} \times \frac{100.00 \text{ mL pickle water}}{1 \text{ pickle}} = 2.2 \times 10^2 \text{ mg Na / pickle}$

$$\frac{44. \text{ mg Na}}{20.00 \text{ mL pickle water}} \times \frac{100.00 \text{ mL pickle water}}{1 \text{ pickle}} = 2.2 \times 10^2 \text{ mg Na / pickle}$$

$$\frac{44. \text{ mg Na}}{20.00 \text{ mL pickle water}} \times \frac{100.00 \text{ mL pickle water}}{1 \text{ pickle}} = 2.2 \times 10^2 \text{ mg Na / pickle}$$

Average:  $2.2 \times 10^2 \text{ mg Na / pickle}$

Standard Deviation:  $0.0 \text{ mg Na / pickle}$

4. The weight of my pickle spear was 32.5753 g. I used a Mt Olive pickle spear, which has 260 mg of sodium. The 95% confidence interval for my data is limited to just  $2.2 \times 10^2 - 2.2 \times 10^2 \text{ mg}$  since the precision of the measurements was lower than their deviation, so twice the standard deviation is still 0.

# Lab Notebook

Exp. No. 4	Experiment/Subject: Mohr Titration for Analysis of Sodium in Pickles	Date: 2/13/25	Course & Section No. 2322
Name: Nathaniel White	Lab Partner: Lukas Rospko	Locker/Desk No.	

  

**Objective:** Determine concentration of Sodium in pickle, compare with listed nutrition facts.

**Procedure:**

**A** Blend a pickle into a liquid, having blotted it dry with a paper towel beforehand. Place this liquid to 100mL with water.

**B** Weigh three 0.1g samples. Place 0.1000g NaCl into 125mL erlenmeyer flasks three times. For each:

- Add 50mL DI & stirbar; mix
- Add sodium bicarbonate slowly until no fizzling.
- Add 2mL of 5% K<sub>2</sub>CrO<sub>4</sub>
- Fill buret with 0.1000M silver nitrate and titrate until red

**C** Take three 20mL aliquots of the pickle water from A and put them into 125mL erlenmeyer flasks. Follow steps in B for each.

Pickle nutritional label  
Mt Olive bottle 280mg Na

Mass pickle: 32.5753 g

W/ OXALATES PRESENT

Wash samples:

	Vol. DI	Vol. K <sub>2</sub> CrO <sub>4</sub>
Std A: 0.1034	50.08mL	2.01mL
Std B: 0.0992	49.71mL	2.03mL
Std C: 0.0982	49.24mL	1.90mL

Volume

Titration A: 24.30mL - 49.24mL  
B: 9.18mL - 26.15mL  
C: 26.39mL - 43.52mL

Pickle samples

	Vol DI	Vol K <sub>2</sub> CrO <sub>4</sub>
1	50.0mL	1.91mL
2	49.2mL	2.00mL
3	49.0mL	2.00mL

Pickle titrations

4: 11.20mL - 30.37mL  
2: 30.68mL - 49.00  
added to 19.30  
- 20.58  
3: 20.84 - 40.27

  

Signature: Nathaniel White	Date:	Witness/TA: [Signature]	Date: 2/15/25
----------------------------	-------	-------------------------	---------------

THE HAYDEN-McNEIL STUDENT LAB NOTEBOOK

Note: Place fold-over back cover under copy sheet before writing