



## Lab 9: Analysis of Wine Using PCA

### Questions

1. [https://docs.google.com/spreadsheets/d/11nICFs\\_IWLvptGVQwJOkuhFIGhMHvmXj-roou\\_IRF8I/edit?usp=sharing](https://docs.google.com/spreadsheets/d/11nICFs_IWLvptGVQwJOkuhFIGhMHvmXj-roou_IRF8I/edit?usp=sharing)
- 2.





### 3. Voltages given relative to SHE

Wine type and name	Scan Rate (mV/s)	E <sub>pa</sub> (mV)	E <sub>pc</sub> (mV)	E <sub>p/2</sub> (mV)	E <sub>mid</sub> (mV)	(E <sub>pa</sub> +E <sub>p/2</sub> )/2 (mV)	E <sub>pa</sub> -E <sub>p/2</sub>	I <sub>pa</sub> (A)	I <sub>pc</sub> (A)	I <sub>pc</sub> /I <sub>pa</sub> (A)
White Italy Pinot Grigio B	20	1419.7	910.5	1297	1165.1	1358.35	122.7	3.25E-06	-3.00E-06	-0.92
Rose Australia	20	1453	913	1308	1183	1380.5	145	3.12E-06	-2.96E-06	-0.95
Unknown	20	1421.6	927	1313	1174.3	1367.3	108.6	-3.06E-06	-3.06E-06	1.00

4. Given that we only have three data points for each category, there is insufficient power to determine based on our samples alone which characteristics are the same or different between the wine samples. All observed and computed values are relatively similar to roughly the same extent. The unknown wine is closer to the white wine for cathodic measurements and closer to the Rose for anodic measurements.

The potential for ferricyanide vs Ag/AgCl can be computed from literature as  $360 \text{ mV} - 220 \text{ mV} = 140 \text{ mV}$ . It was not recorded experimentally during this lab due to explicit instruction by a TA.

Our blank reading showed strong repeated evidence of contamination with observed peaks similar to our voltammograms for the wine sample, thus we did not subtract our blank for our data analysis.

[https://docs.google.com/spreadsheets/d/10stXnoedHNNH2xuzVSs6tJZiiTX57NDyTO8v6CKn\\_eck/edit?gid=301451791#gid=301451791](https://docs.google.com/spreadsheets/d/10stXnoedHNNH2xuzVSs6tJZiiTX57NDyTO8v6CKn_eck/edit?gid=301451791#gid=301451791)

# Lab Notebook

Exp. No. 9	Experiment/Subject Analysis of wine using cyclic voltammetry and PCA	Date Apr 3
Name Neil Hines / Wynn	Lab Partner Lalces Rosendo	Locker/Desk No.
		Course & Section No.

Safety: Ethanol is very flammable, strong acids and bases are corrosive

Objective: Get CV scans for wines to determine identity of unknown with PCA

Procedure

A: <sup>Get</sup> make 100 mL  $\text{KNO}_3$  @ 1M

~~use 10.1032 g~~

Make 10 mL  $\text{K}_3\text{Fe}(\text{CN})_6$

soln with 0.0330 g

$\text{K}_3\text{Fe}(\text{CN})_6$  in 1M  $\text{KNO}_3$  soln

→ 0.0353 g  $\text{K}_3\text{Fe}(\text{CN})_6$

C

White Wine: Barbaresco Italy Pinot Grigio

Rosé: Australia

Both were diluted  $\frac{100 \mu\text{L}}{10 \text{ mL}}$

D CV

Platinum electrode,  $\text{K}_3\text{Fe}(\text{CN})_6$

Scan rate 100 mV/s  $\text{K}_3\text{Fe}(\text{CN})_6$

100 mV - 400 mV

Build electrode, 300 - 700 mV

4 seg meters, discard first 2

20 mV/s scan rate

used for wines

Signature	Date Apr 3	Witness/TA	Date 4/3/25
-----------	------------	------------	-------------

THE HAYDEN-McNEIL STUDENT LAB NOTEBOOK

Note: Place fold-over back cover under copy sheet before writing