



## Lab 9: Analysis of Wine Using PCA

### Questions

1. [https://docs.google.com/spreadsheets/d/11nICFs\\_IWLvptGVQwJOkuhFIGhMHvmXj-roou\\_IRF8I/edit?usp=sharing](https://docs.google.com/spreadsheets/d/11nICFs_IWLvptGVQwJOkuhFIGhMHvmXj-roou_IRF8I/edit?usp=sharing)

2.





### 3. Voltages given relative to SHE

Wine type and name	Scan Rate (mV/s)	E <sub>pa</sub> (mV)	E <sub>pc</sub> (mV)	E <sub>p/2</sub> (mV)	E <sub>mid</sub> (mV)	(E <sub>pa</sub> +E <sub>p/2</sub> )/2 (mV)	E <sub>pa</sub> -E <sub>p/2</sub>	I <sub>pa</sub> (A)	I <sub>pc</sub> (A)	I <sub>pc</sub> /I <sub>pa</sub> (A)
White Italy Pinot Grigio B	20	1419.7	910.5	1297	1165.1	1358.35	122.7	3.25E-06	-3.00E-06	-0.92
Rose Australia	20	1453	913	1308	1183	1380.5	145	3.12E-06	-2.96E-06	-0.95
Unknown	20	1421.6	927	1313	1174.3	1367.3	108.6	-3.06E-06	-3.06E-06	1.00




4. Given that we only have three data points for each category, there is insufficient power to determine based on our samples alone which characteristics are the same or different between the wine samples. All observed and computed values are relatively similar to roughly the same extent. The unknown wine is closer to the white wine for cathodic measurements and closer to the Rose for anodic measurements.

The potential for ferricyanide vs Ag/AgCl can be computed from literature as  $360 \text{ mV} - 220 \text{ mV} = 140 \text{ mV}$ . It was not recorded experimentally during this lab due to explicit instruction by a TA.




Our blank reading showed strong repeated evidence of contamination with observed peaks similar to our voltammograms for the wine sample, thus we did not subtract our blank for our data analysis.

## Conclusion




- Regardless of the choice of all features, just the reduction features, or just the oxidation features, the clustering cleanly appears around the lab group which conducted the analysis, not the color or country of the wine. These clusters have been circled and labeled with the initials of the lab group which contributed the sample.



**OVERALL**



# OXIDATION



## REDUCTION

- Our results are inconclusive. Clustering based on partner pair rather than the actual wine sample indicates the noise from instrumentation or operator error greatly outweighs the signal. We don't see in any of the observation plots the ability to neatly group data points based on color or country of origin. We also don't see the expected result of Rose between in between white wines and red wines; rather, we see that the same Rose analyzed by three different groups does not appear grouped together in any plot.
- Both analyses were inconclusive, so no observations can be made. However, univariate analysis is generally more sensitive when there is an expected pattern, whereas multivariate analysis can identify unforeseen patterns.

# Lab Notebook

Exp. No. 9	Experiment/Subject Analysis of wine using cyclic voltammetry and PCA	Date Apr 3
Name Neil Hines/Whitney	Lab Partner Talia Rose	Locker/Desk No.
		Course & Section No.

Safety: Ethanol is very flammable, strong acids and bases are corrosive

Objective: Get CV scans for wines to determine identity of unknown with PCA

Procedure

A: <sup>Get</sup> make 100 mL  $\text{KNO}_3$  @ 1M

~~use 10.1032 g~~

Make 10 mL  $\text{K}_3\text{Fe}(\text{CN})_6$

soln with 0.0330 g

$\text{K}_3\text{Fe}(\text{CN})_6$  in 1M  $\text{KNO}_3$  soln

→ 0.0353 g  $\text{K}_3\text{Fe}(\text{CN})_6$

C

White Wine: Barbera Italy Pinot Grigio

Rosé: Australia

Both were diluted  $\frac{100 \mu\text{L}}{10 \text{ mL}}$

D CV

Platinum electrode,  $\text{K}_3\text{Fe}(\text{CN})_6$

Scan rate 100 mV/s  $\text{K}_3\text{Fe}(\text{CN})_6$

100 mV - 400 mV

Build electrode, 300 - 1400 mV

4 seg meters, discard first 2

20 mV/s scan rate

used for wines

Signature	Date Apr 3	Witness/TA	Date 4/3/25
-----------	------------	------------	-------------

THE HAYDEN-McNEIL STUDENT LAB NOTEBOOK

Note: Place fold-over back cover under copy sheet before writing