

## Lab 7: Analysis of Caffeine in an OTC Analgesic by HPLC

### Questions

1. Performing linear regression on the peak area from our caffeine standard calibration curve yields a slope of 11169.7 AU/ $\mu$ L stock solution and a y-intercept of 70876.1 AU. Our area for the peak corresponding to caffeine for anacin had an area of 3228713 AU, which corresponds to 282.714  $\mu$ L stock solution on the calibration curve. Some simple dimensional analysis yields that  
$$\frac{282.714 \mu\text{L stock solution}}{1000 \mu\text{L anacin solution}} * 250 \text{ mL anacin solution} = \frac{1}{4} * 282.714 \text{ mL stock solution} = \frac{1}{4} * 282.714 \text{ mL stock solution} * \frac{12.4 \text{ mg caffeine}}{25 \text{ mL stock solution}} = 35.1 \text{ mg caffeine in the anacin tablet.}$$
2. The sample was 507.6 mg, so the tablet is 6.91% caffeine w/w.
3. N/A, only one replicate
4. Result seems reasonably close to the expected value of 32 mg.

# Lab Notebook

Exp. No. <u>7</u>	Experiment/Subject <u>Analysis of an OTC Analgesic by HPLC</u>	Locker/ Desk No.	Course & Section No.
Name <u>Nathanial White</u>	Lab Partner <u>Lukas Rosicko</u>		

  

**Objective:** Use HPLC to quantify aspirin and caffeine in an Analgesic Tablet

**Procedure:**

Caffeine stock solution  
 stock soln. 12.5mg caffeine in 25ml v.f.  
 1ml with HPLC water

Aspirin standards, with 50  $\mu$ L, 150  $\mu$ L, 200  $\mu$ L, and 300  $\mu$ L of stock, rest HPLC water

Analgesic Mixture

- one Analgesic tablet was ground to a fine powder
- dissolved in HPLC water with sonicator
- After a large amount had dissolved, 250 mL flask filled with HPLC water
- Transferred through nylon syringe filter before being aliquotted into autosampler vials

Aspirin standard: 1.5 mL @ 0.0010 g/mL (premade)

Other PR strips - 0.0124 g Caffeine Flinn Lot 92331

used Fisher Brand #18 for micropipetting

Acculab Analytic Balance S/N 11904121  
 - 0.5076 g Analgesic

HPLC

- A8 - Blank
- B1 - 50
- B2 - 150
- B3 - 200
- B4 - 300
- B5 - Aspirin
- B6 - analgesic

  

Signature <u>Nathanial White</u>	Date	Witness TA <u>Lukas Rosicko</u>	Date
----------------------------------	------	---------------------------------	------

THE HAYDEN-McNEIL STUDENT LAB NOTEBOOK Note: Place fold-over back cover under copy sheet before writing